

Research Core Unit Transcriptomics

Crude probe characterization workflow

Version: 22.08.2012

Introduction

The crude probe characterization workflow is performed by us to allow for a rapid discrimination between probes that solely serve as technical controls and probes, directed against endogenous transcripts. Furthermore, the degree of functional characterization and the annotation state of the respective transcripts (genes) is considered and used to assign a specific attribute to all probes, according to the following 3 possibilities* :

- “1) gene symbol”
- “2) poorly characterized”
- “3) control”

The source of annotation data for Agilent’s microarray probes

Annotation data is provided at Agilent’s eArray portal (Figure 1). Regularly updated annotation files (Figure 2) are downloaded by us from this site and serve as basis for the annotation block within our standard Excel result files and for the crude probe characterization workflow (as described below).

* Please also note paragraph „ Important explanatory information “ at the end of this manual.



Figure 1: The eArray login page.

1	Column	Row	Name	ID	RefNumber	ControlType	GeneName	TopHit	Description	Go	Chromosoma	EntrezGeneID
2		1	GE_BrightCo	GE_BrightCo	1	pos	GE_BrightCorner		Unknown		Unknown	
3		2		DarkCorner	2	pos	DarkCorner		Unknown		Unknown	
4		3		DarkCorner	3	pos	DarkCorner		Unknown		Unknown	
5		4		DarkCorner	4	pos	DarkCorner		Unknown		Unknown	
6		5		DarkCorner	5	pos	DarkCorner		Unknown		Unknown	
7		6		DarkCorner	6	pos	DarkCorner		Unknown		Unknown	
8		7		DarkCorner	7	pos	DarkCorner		Unknown		Unknown	
9		8		DarkCorner	8	pos	DarkCorner		Unknown		Unknown	
10		9		DarkCorner	9	pos	DarkCorner		Unknown		Unknown	
11		10		DarkCorner	10	pos	DarkCorner		Unknown		Unknown	
12		11		DarkCorner	11	pos	DarkCorner		Unknown		Unknown	
13		12	1	NM_152565	A_23_P1461	12	false	ATP6V0D2	ref NM_1525 Homo sapien GO:0005515 chr8:871653 entg 245972			
14		13	1	NM_004333	A_23_P4293	13	false	BRAF	ref NM_0043 Homo sapien GO:0000166 chr7:140434 entg 673			
15		14	1	NM_015987	A_23_P1170	14	false	HFBP1	ref NM_0159 Homo sapien GO:0005488 chr12:13127 entg 50865			

Figure 2: An example of an Agilent annotation file (AMADID 026652) available at the eArray portal. Such files serve as basis for the annotation block within our standard Excel result files and for the crude probe characterization workflow (Table 1).

As depicted in Figure 3, one of three possible attributes is assigned to each probe and is integrated into our standard Excel result file within column: “crude probe characterization (1-3)” (1-3)”.

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE	AF	AG	AH	AI			
1	S	S	S	S	S	S	S	S	S	S	S	S	S	S	S	S	S	S	Fea	RS	crude probe characterization (1-3)	num	GeneName	D	gls	M3622	c	gls	M3623	c	gls	M3624	c	gls	M3625	c	gls	
2	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	1	(GE	(3)	control	14	GE_BrightCorner	U	1	51788	0	1	74496	0	1	81205	0	1	76466	0	
3	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	2	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
4	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	3	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
5	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	4	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
6	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	5	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
7	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	6	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
8	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	7	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
9	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	8	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
10	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	9	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
11	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	10	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
12	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	11	DarkC	(3)	control	604	DarkCorner	U	0	15	0	0	15	0	0	15	0	0	15	0	0
13	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	12	AN	re	(1)	gene symbol	1	ATP6V0D2	H	1	47	0	1	43	0	1	34	0	1	37	0
14	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	13	AN	re	(1)	gene symbol	10	BRAF	H	1	315	0	1	315	0	1	297	0	1	332	0
15	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	14	AN	re	(1)	gene symbol	10	HEBP1	H	1	7963	0	1	7426	0	1	6652	0	1	8588	0
16	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	15	AN	re	(1)	gene symbol	10	RPAP3	H	1	2517	0	1	2532	0	1	2283	0	1	2712	0
17	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	16	A	re	(2)	poorly characterized	1	A_24_P358131	U	1	799	0	1	717	0	1	676	0	1	693	0

Figure 3: A part of the standard result file for single-color studies. The “crude probe characterization (1-3)” and the “GeneName” columns are highlighted in red.

Crude probe characterization workflows

The following workflow is used to assign the described attributes to all probes of Agilent's mRNA expression microarrays utilized by us. The column names refer to the annotation file as exemplified in Figure 2.

Table 1: The crude probe characterization workflow for human microarrays of design types AMADID 014850 and 026652, for murine microarrays of design types AMADID 014868 and 026655, and for rat microarrays of design types AMADID 014879 and 028282.

Step	Column to check in annotation file	Search pattern	Probe characterization
1		set all probes to	2) poorly characterized
2	Name	"NM_*"	1) gene symbol
3	Description	"*hypothetical*" or "*predicted*" (case insensitive), "*RIKEN*" or "*FLJ#*" or "*Mus musculus expressed sequence*" or "*Mus musculus cDNA sequence*"	2) poorly characterized
4	GeneName	"KIAA*" or "LOC*" or "hCG_*" or "RP11-*" or "RP1 *" or "RP3-*" or "RP6-*" or "MGC*" or "RGD*", "CXorf#*" or "CX#orf#*" or "CYorf#*" or "CY#orf#*" or "C#orf#*" or "C##orf#*" or "CXORF#*" or "CX#ORF#*" or "CYORF#*" or "CY#ORF#*" or "C#ORF#*" or "C##ORF#"	2) poorly characterized
5	ControlType	not "false"	3) control
6		if one replicate probe of a gene is characterized as "1) gene symbol" then all other replicates are also set to	1) gene symbol

- exactly one number is expected (0-9)

* - zero or more letters, numbers or any other character are allowed

Important explanatory information

Please note, that the discrimination between "1) gene symbol" and "2) poorly characterized" in its present form is not completely conclusive in a semantical sense. According to our workflow, we decided to allocate the attribute "1) gene symbol" only for protein-coding transcripts (genes). Thus, all non-coding transcripts (e.g. those with an "NR_*" RefSeq accession entry) receive the attribute "2) poorly characterized", irrespective of their actual characterization/annotation status (see Figure 4 for an example).

On the other hand, there are some examples for transcripts (genes) for which no RefSeq accession entry of the type: "NM_*", but only an ENSEMBLE database entry exists: accession entry of the type "ENST*" (human), "ENSMUST*" (mouse) or "ENSRNOT*" (rat). Anyhow, even some of these transcripts (genes) could have a gene symbol officially allocated (Figure 5). In such a situation we decided to assign the respective probes with the attribute "2) poorly characterized", unless, additional probes, directed against the same transcript (gene) and possessing an "NM_*" accession entry are present on the microarray (see Table 1, steps 2 and 6).

